
Liquid Needle: revolutionary dosing technique for contact angle scientifically validated

- Scientific study confirms accuracy of results compared with classic needle dosing
- User-independent dosing dynamic reliably prevents unwanted pre-wetting

With the Liquid Needle dosing technique, KRÜSS GmbH is since recently using a new method of drop positioning for measuring static contact angle. The pressure-based method significantly accelerates the dosing process compared with classic needle dosing. Here, dosing units arranged in parallel produce two drops at the same time enabling measurements of the surface free energy to be carried out within one second. A current scientific study in the journal *Colloid and Polymer Science* verifies that contact angles are measured accurately in spite of this high speed.

During the study, contact angles were measured on 14 different materials with the new pressure dosing method and the classic needle dosing technique. With a wide range of hydrophobic and hydrophilic, rough and smooth as well as chemically pure and technical surfaces, the investigation covers all practical cases. The results of the scientifically reviewed publication show without exception a good accord between the contact angle results of the two dosing methods.

Alternatives to needle dosing have previously failed due to the region outside the contact area of the drop being pre-wetted as the dynamics were too high. This falsified the contact angle resulting in values which were too small. In contrast, with the KRÜSS Liquid Needle method, drops are formed in a controlled manner and with low dynamics with the help of a continuous jet. The jet is so thin compared with the final drop size that the contact area does not spread to any greater extent than with needle dosing. This is clearly shown by the results of the study.

Both techniques are therefore basically distinguished by careful dosing. However, this is not always the case with classic needle dosing. The user can inadvertently increase the dynamics and thus falsify the value by an unsuitable choice of depositing speed or dosing distance. This is ruled out technically with the new dosing method. The Liquid Needle is also particularly beneficial with liquid-repellent samples. While smaller drops with needle dosing are difficult to deposit on such hard-to-wet materials, dosing with the Liquid Needle presents no difficulties.

Overall, the study shows that the Liquid Needle combines speed with high accuracy, reliably prevents unwanted pre-wetting, and furthermore enormously simplifies dosing on liquid-repellent samples. So far, KRÜSS is using the new dosing technique in the mobile test instrument *Mobile Surface Analyzer – MSA* and making it available for all stationary contact angle measuring instruments.

Literature:

Ming Jin, Raymond Sanedrin, Daniel Frese, Carsten Scheithauer, Thomas Willers: Replacing the solid needle by a liquid one when measuring static and advancing contact angles. Colloid and Polymer Science, DOI 10.1007/s00396-015-3823-1 (2016).

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Literatur: *Ming Jin, Raymond Sanedrin, Daniel Frese, Carsten Scheithauer, Thomas Willers: Replacing the solid needle by a liquid one when measuring static and advancing contact angles. Colloid and Polymer Science, DOI 10.1007/s00396-015-3823-1 (2016).* ■

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31ST INTERNATIONAL EXHIBITION FOR FINE AND SPECIALITY CHEMICALS

The fine & speciality chemicals exhibition

Chemspec Europe is exclusively dedicated to the fine, custom and speciality chemicals sector. With its focussed profile it is the place to be for international industry professionals. Leading manufacturers, suppliers and distributors will present a wide spectrum of fine and speciality chemicals for various industry sectors.

Establish new business contacts, benefit from excellent networking opportunities and be inspired by the latest results in Research and Development at top-class conferences. Be part of the show where ideas become innovations!

TOP-CLASS CONFERENCES & WORKSHOPS!

- Royal Society of Chemistry Symposium
- REACHReady Regulatory Services Conference
- Agrochemical Intermediates Conference
- The Pharma Outsourcing Best Practices Panel

www.chemspeceurope.com

1 - 2 JUNE 2016
BASEL MESSE, SWITZERLAND

Organisers: Mack Brooks Exhibitions Ltd